
# LIWO HemaStainer

Click & Walk Away

Stained, Washed & Ready in 90 Seconds






With Livo HemaStainer, gone are the days when your only option for processing blood samples was relying in specialised manual labor, a luxury in terms of turnaround times that not all labs and healthcare practices can afford.

Designed with the perishable nature of blood products in mind, this compact, fully automated device measures only 31.5 x 28 x 18cm and allows you to stain any number of blood smears with unrivalled speed and precision, eliminating idle time and the risk of cross contamination.

#### "Two-step, single-click operation."

**Step 1:** Load your smeared slide into the Livo HemaStainer.

**Step 2:** Push the button. The device dispenses a precise amount of fresh stain, stains the blood sample, and washes the slide, all in exactly 90 seconds. Once the process is finished, the HemaStainer emits a beep so that you can take the slide for analysis and load a new one.



## In Pursuit of Highest Staining Quality

#### Textbook quality staining at the push of a button

The automated two-step process ensures consistent staining quality every time, without fail.

#### No cross contamination

Livo HemaStainer uses fresh stain every time, thus eliminating the risk of cross contamination.

### Cost Savings

Maximum efficiency & minimal stain usage



Apart from being the most affordable automated hematology stainer on the market, the Livo HemaStainer allows you to save thousands of rupees in yearly material costs due to its precise use of stain for each slide, resulting in zero waste.

#### **Save Time**

Stain in 90 Seconds & shorten TAT



Unlike the typical 10-12 minute staining time, the Livo HemaStainer stains and washes your slides in a mere 90 seconds with just the push of a button, reducing both process steps and idle time common to manual staining.

